



CiDehol[®]ST₇₀

Technical Data Sheet

-BRIEF DESCRIPTION -

CiDehol ST is a sterile 70% (v/v) USP isopropyl alcohol solution made with Water for Injection (WFI). It can be used to clean surfaces and equipment found in cleanroom areas within pharmaceutical, biotechnology, and medical device facilities or anywhere the use of a sterile, non-pyrogenic alcohol solution is required.

CiDehol ST is filtered to 0.22 microns, bottled and double-bagged within a CLASS 100 cleanroom. It is gamma-irradiated to a sterility assurance level (SAL) of 10⁻⁶. Each lot is tested for sterility via USP 14-day sterility test and for bacterial endotoxin via the LAL test prior to release. A copy of lot specific documentation, which details QC, irradiation, sterility, and bacterial endotoxin limits, is attached to every case of CiDehol ST. All lots of Cidehol ST are guaranteed to contain less than 0.25 EU/mL, the endotoxin level criterion for WFI.

- USE PROTOCOL -

Use in well-ventilated area. Hold spray bottle six to eight inches from surface and spray area until it is completely covered with solution. Or wipe CiDehol ST onto surface using cleanroom grade wipers. Allow to air dry.

- INGREDIENTS -

70% Isopropyl Alcohol, USP
30% Water for Injection, USP
by Volume

- SHELF LIFE AND STORAGE -

The shelf life is 2 years. The expiration date is on the product. Do not store at temperatures above 120°F.

- PRECAUTIONARY STATEMENT -

Flammable - do not come in contact with flame or extreme heat. Avoid eye and skin contact and ingestion.

- PACKAGING -

16 oz. Trigger Spray / 12 per case
1 Gallon / 4 per case

- CATALOG NUMBERS -

<u>Manufacturer</u>	<u>Item</u>	<u>Cat. No.</u>
Decon Labs, Inc.	16 oz. Bottle	8316
	1 Gal Bottle	8301

**For Additional Technical Information call
Decon Customer Service
(800) 332-6647, or
See our website at www.deconlabs.com**

®Decon Labs, Inc.

Revised 9/21/2011



CiDehol[®]ST₇₀

Technical Data Sheet

-BRIEF DESCRIPTION -

CiDehol ST is a sterile 70% (v/v) USP isopropyl alcohol solution made to clean surfaces and equipment found in cleanroom areas within pharmaceutical, biotechnology, and medical device facilities or anywhere the use of a sterile, non-pyrogenic alcohol solution is required.

CiDehol ST is filtered to 0.22 microns, bottled and double-bagged within a CLASS 100 cleanroom. It is gamma-irradiated to a sterility assurance level (SAL) of 10⁻⁶. Each lot is tested for sterility via USP 14-day sterility test and for bacterial endotoxin via the LAL test prior to release. A copy of lot specific documentation, which details QC, irradiation, sterility, and bacterial endotoxin limits, is attached to every case of CiDehol ST. All lots of Cidehol ST are guaranteed to contain less than 0.25 EU/mL, the endotoxin level criterion for WFI.

- USE PROTOCOL -

Use in well-ventilated area. Hold spray bottle six to eight inches from surface and spray area until it is completely covered with solution. Or wipe CiDehol ST onto surface using cleanroom grade wipers. Allow to air dry.

- INGREDIENTS -

70% Isopropyl Alcohol, USP
30% Purified Water, USP
by Volume

- SHELF LIFE AND STORAGE -

The shelf life is 2 years. The expiration date is on the product. Do not store at temperatures above 120°F.

- PRECAUTIONARY STATEMENT -

Flammable - do not come in contact with flame or extreme heat. Avoid eye and skin contact and ingestion.

- PACKAGING -

12 oz. Trigger Spray / 16 per case
32 oz/ 12 per case

- CATALOG NUMBERS -

<u>Manufacturer</u>	<u>Item</u>	<u>Cat. No.</u>
Decon Labs, Inc.	12oz. Bottle	8312
	32oz Bottle	8332

For Additional Technical Information call
Decon Customer Service
(800) 332-6647, or
See our website at www.deconlabs.com

®Decon Labs, Inc.

Revised 9/21/2011



CiDehol® ST 70

Hands-Free Dispenser

Technical Data Sheet

- BRIEF DESCRIPTION -

CiDehol ST is a sterile 70% (v/v) USP isopropyl alcohol solution made to clean gloved hands in gowning areas and cleanroom areas within pharmaceutical, biotechnology, and medical device facilities or anywhere the use of a sterile, non-pyrogenic alcohol solution is required.

CiDehol ST is filtered to 0.22 microns, filled and double-bagged within a CLASS 100 cleanroom. It is gamma-irradiated to a sterility assurance level (SAL) of 10⁻⁶. Each lot is tested for sterility via USP 14-day sterility test and for bacterial endotoxin via the LAL test prior to release. A copy of lot specific documentation, which details QC, irradiation, sterility, and bacterial endotoxin limits is attached to every case of CiDehol ST. The designation WFIQ indicates that CiDehol ST is made with USP grade Purified Water which undergoes bacterial endotoxin testing. All lots of Cidehol ST are guaranteed to contain less than 0.25 EU/mL, the endotoxin level criterion for WFI.

- USE PROTOCOL -

Decon 8329 in 32 oz. bags are to be used with Touch-free Dispenser Part Number 8308. Touch Free dispenser operates on infrared technology. The dispenser is mounted via hardware or adhesive strips supplied with the unit. It can be battery (four size “D” alkaline batteries-not included with unit) or A/C operated, power supply is included. Follow the user manual for installing the bags.

- INGREDIENTS -

70% Isopropyl Alcohol, USP
30% Purified Water, USP
by Volume

- SHELF LIFE AND STORAGE -

The shelf life is 2 years. The expiration date is on the product. Do not store at temperatures above 120°F.

- PRECAUTIONARY STATEMENT -

Flammable - do not come in contact with flame or extreme heat. Avoid eye and skin contact and ingestion.

- PACKAGING -

32 oz. Sterile Bag / 8 per case*

***For use with Dispenser Decon Part # 8308**

- CATALOG NUMBERS -

<u>Manufacturer</u>	<u>Item</u>	<u>Cat. No.</u>
Decon Labs, Inc.	32 oz. (Bag)	8329
	Dispenser	8308

**For Additional Technical Information call
Decon Customer Service
(800) 332-6647, or
See our website at www.deconlabs.com**

®Decon Labs, Inc.

Revised 09/21/2011