



ELIMINase®

RNase/DNase Remover



ELIMINase® removes RNase, DNase, unwanted RNA, unwanted DNA and radioactive contamination. It works instantly so there is no waiting. ELIMINase is an ideal substitute for using DiEthylPryoCarbonate (DEPC) which is a carcinogen. It also eliminates the time-consuming method of autoclaving. **ELIMINase®** can be used on all surfaces including glass, plastic, metal, and ceramic. ELIMINase is phosphate-free, safely disposable, non-toxic, and totally rinsable for the decontaminated surface.

USE PROCEDURES: Use **ELIMINase®** full strength.

METHOD 1

Thoroughly wet item to be decontaminated with **ELIMINase®** - ensure all surfaces are coated. Swirl for a few seconds, drain, and immediately rinse twice with the RNase-free water.

METHOD 2

Soak item in **ELIMINase®** until needed. Drain and rinse as above.

METHOD 3

Soak an RNase-free wipe in **ELIMINase®**, wipe over surface to be decontaminated. Rinse as above.

For total radioactive decontamination, soak item in **ELIMINase®** for 2-12 hours. Rinse as above. Wipe over benchtops and sinks to remove residual contamination.

CAUTION:

Alkaline! Severe eye irritant, mild skin irritant. Wear eye and skin protection when handling **ELIMINase®**.

Packaging

Size	Case QTY	Decon #
250mL	12	1101
950mL	12	1102

Rev. 12/2022