



Hyde-Away™

NEUTRALIZING LIQUID



HydeAway™ Neutralizing Liquid is a very unique product that neutralizes and crosslinks used aldehyde into a non-toxic, non-hazardous end product in a 1:15 ratio. The 1 gallon container of HydeAway™ contains sufficient liquid to convert a total of fifteen (15) gallons in either single or multiple dose applications. Treating used aldehyde with HydeAway™, prior to drain disposal, is better for the environment and makes for a safer workplace.

HydeAway™ makes regulatory compliance easy, inexpensive and convenient.

USE PROTOCOL: Wear hand, eye and respiratory protection. Transfer aldehyde waste to appropriate disposal container. Slowly pour in the HydeAway™. The transfer should be performed in a well ventilated area or under a hood. The neutralizing process produces a gas that will visibly “bubble” as the gas is released. The container should remain open until the neutralization is complete. When used as directed, eight (8) ounces, or two hundred fifty (250) mL. of

HydeAway™ will treat one gallon of used aldehyde. Allow a minimum of one hour for completion of the neutralization treatment process prior to final disposal. General laboratory safety practices must be followed at all times. When used as directed, HydeAway™ converts aldehyde to a nontoxic and non-hazardous material.

STORAGE & DISPOSOL: Safe storage and disposal of chemicals is an integral part of any laboratory operation. Store in original container with lid tightly sealed. Keep from heat and open flame. Properly dispose of aldehyde treated with HydeAway™ according to institutional procedures, as well as local, state and federal regulations. Many areas allow used aldehyde treated with HydeAway™ to be placed directly into the drain, the regular solid waste system or incinerated.

PACKAGING

4 X 1 GALLON

Decon Catalog No. 2201

Rev. 09/2022