



HydroSolv™

Solvent and Clearing Agent

Technical Data Sheet

- BRIEF DESCRIPTION -

HydroSolv is a clearing agent and dehydrating solvent to be used in histological procedures. It is a safe substitute for xylene for the production of the highest quality tissue specimens. HydroSolv is a unique isoparaffinic solvent blend with a flash point that makes it nonflammable (combustible). HydroSolv will not remove bound water and hence it will keep tissues in a more natural state. HydroSolv is non-irritating, non-sensitizing, greaseless and practically odorless. HydroSolv contains no limonene or citrus oils.

- USE PROTOCOL -

Use for manual or automatic tissue processor staining procedures in histology/cytology applications that require xylene.

- STORAGE AND DISPOSAL -

Warning. HydroSolv is combustible. Keep from heat and open flame. Store in original container with lid tightly sealed. Follow federal, state and local restrictions for disposal.

- PRECAUTIONARY STATEMENT -

Keep out of reach of children.
Avoid contact with eyes, skin, and clothing. Avoid ingestion and inhalation. Empty containers retain product residue and can be

dangerous. Avoid contact with heat, sparks, or flame.

FIRST AID:

Eye Contact: Immediately flush eyes with plenty of water for at least fifteen minutes, occasionally lifting the upper and lower eyelids. Contact a physician.

Skin Contact: Flush skin with plenty of soap and water for at least fifteen minutes while removing contaminated clothing and shoes. Contact a physician if irritation develops or persists.

Ingestion: Do NOT induce vomiting. Never give anything by mouth to an unconscious person. Contact a physician immediately.

Inhalation: Remove from exposure to fresh air immediately. Contact a physician.

- PACKAGING -

1 Gallon / 4 per case

- CATALOG NUMBERS -

<u>Manufacturer</u>	<u>Item</u>	<u>Cat. No.</u>
Decon Labs, Inc.	1 Gal.	2301

For Additional Technical Information call
Decon Customer Service
(800) 332-6647, or
See our website at www.deconlabs.com

™Decon Labs, Inc.

Revised 11/16/06