



JET STREAM™

High Pressure Dust Remover

Technical Data Sheet

- BRIEF DESCRIPTION -

Jet Stream is moisture-free propelled air used to remove microscopic dust from negatives, microscope slides, lenses, balances and other laboratory equipment.

use of product, or leave the area. Avoid breathing vapors, avoid contact with skin and eyes.

- USE PROTOCOL -

Remove tamper-resistant tab on top of trigger sprayer. Insert tube firmly into front of sprayer. Hold can 2-4 inches from surface while operating. Do not use this product on camera mirrors. Always use this product in an upright position and do not shake it. Do not let product in liquid form come in contact with film negatives or other samples.

First Aid: In case of eye contact, flush thoroughly with plenty of water for fifteen minutes and get medical attention. For liquid contact on skin, check for frostbite, contact a physician if burn occurs. If inhaled and symptoms occur, remove to fresh air and get medical attention.

INTENTIONAL MISUSE OF THIS PRODUCT BY DELIBERATELY CONCENTRATING AND INHALING THE CONTENTS MAY BE HARMFUL OR FATAL!

- STORAGE AND DISPOSAL -

Exposure to temperatures above 120°F may cause can to burst with violence and cause injury. Therefore, do not leave can in enclosed vehicles, direct sunlight, hot water, or near any heat sources. Do not puncture or incinerate, or place in a trash compactor.

- PACKAGING -

11 oz. (311 grams) Aerosol Can
6 per case

- PRECAUTIONARY STATEMENT -

Liquid contents may cause frostbite on contact with skin, therefore do not tilt, shake, or turn can upside down before or during use.

- CATALOG NUMBER -

<u>Manufacturer</u>	<u>Item</u>	<u>Decon Cat. No.</u>
Decon Labs, Inc.	11 oz. (311 grams)	3711

Use in well-ventilated areas. Vapor reduces oxygen available for breathing. If you experience eye watering, headaches, dizziness, increase the fresh air, discontinue

**For Additional Technical Information call
Decon Customer Service
(800) 332-6647, or
See our website at www.deconlabs.com**

®Decon Labs, Inc.

Revised 11/2006