

SPECIALTY CLEANERS



Detergent

Alkaline

Neutral

Acidic

Enzymatic

Concentrated

Ready-to-use

Manual

Soak

Ultrasonic



www.deconlabs.com

DECON'S SPECIALTY CLEANERS FOR ALL CRITICAL APPLICATIONS

SOAK, MANUAL SCRUB, ULTRASONIC CLEANING PH: ALKALINE

CONTRAD® 70 liquid detergent concentrate gets everything clean by soaking, scrubbing and ultrasonics. It rinses away totally, leaves zero residues. Phosphate-free, disposable straight to drain, Contrad 70 is perfect for all critical cleaning, from labware to process equipment, ultrasonic cleaning and radioactive decontamination. The original and still the best... use it and see!



DRI-CLEAN™ is a powdered detergent for soak, manual, and ultrasonic cleaning applications. User friendly bottle provides easy dispensing and eliminates caking. Completely soluble, free-rinsing, and residue-free. Versatile formula removes a wide variety of tenacious contaminants.



NEUTRAD® liquid detergent concentrate for hand scrubbing your labware offers high-performance cleaning with the safety of a near-neutral pH. Phosphate-free, totally rinsable, disposable straight to drain, Neutrad is the perfect replacement for old-fashioned caustic soap powders.



CITRAD™ citric acid-based detergent concentrate is ideal for removing salts, scale, and metal-oxide. Also used for metal brightening and pharmaceutical-process equipment. Completely phosphate-free, totally free-rinsing and biodegradable.



DESCAL® inorganic acid-based detergent and descaler concentrate for removing hard water mineral deposits, rust and crusty build up from labware, water baths and stills. Low-foaming formula allows it to be used for automatic washing applications.



NO-COUNT® surface cleanser is a convenient aerosol foam for removing radioactive contamination from benches, sinks, pipettors, gel boxes, shields and other hard surfaces. Just spray on, then wipe or rinse clean.



ELIMINase® is a non-toxic, ready-to-use solution for instantly removing RNase, DNase, unwanted DNA, and radioactive contamination from labware and equipment. Safe and effective. Eliminates need for autoclaving or using DEPC.



Enzyte™ medical enzyme cleaner is a multi-enzyme detergent concentrate formulated to be especially effective in the removal of proteins, blood, fatty tissues and other organic materials found on medical instruments and labware. Unique triple enzyme formula removes all bioburden and is safe on all instruments and scopes. Non-caustic, non-toxic and disposable straight to drain.



AUTOMATIC WASHERS PH: ALKALINE

CONTRAD® NF non-foaming super-concentrated liquid detergent is formulated for automatic labware washers. Phosphate-free, and totally free-rinsing, Contrad NF maximizes the efficiency of your labware machine.



DRI-CONTRAD® powdered detergent supercharges your automatic labware washer. Non-foaming, instantly soluble, rinsable, and chlorinated for biological control. Dri-Contrad is the only powder you need for your labware washer.



CONTRAD® 100 is a low-foaming detergent designed for clean-in-place systems and other high-pressure spray applications, including automatic washers. Alkaline formulation is phosphate-free, residue-free, and biodegradable.



AVAILABILITY

Item	Size	Decon Cat. #	Case Quantity
Contrad® 70 Soak Cleaner	1 Liter	1002	12 x 1 L
	5 Liters	1003	4 x 5 L
	20 Liters	1004	Each
	210 Liters	1005	Each
Dri-Clean™ Powdered Labware Detergent	2 kG	2502	6 x 2 kG
	50 kG	2501	Each
Neutrad® Manual Cleaner	950 mL	3003	12 x 950 mL
	5 Liters	3001	4 x 5 L
	20 Liters	3002	Each
Citrad® Acid Cleaner & Detergent	3.8 Liters	4401	4 x 3.8 L
DeSCAL® Descaler Concentrate	1 Liter	4001	12 x 1 L
	10 Liters	4010	2 x 10 L
	20 Liters	4002	Each
	115 Liters	4030	Each
No-Count® Radioactive Decontaminant Spray	18 oz.	6518	12 x 18 oz.
ELIMINase® RNase Decontaminant	250 mL	1101	12 x 250 mL
	950 mL	1102	12 x 950 mL
Enzyte™	3.8 Liters	1901	4 x 3.8 L
Contrad® NF Liquid Autowasher Detergent	3.8 Liters	6002	4 x 3.8 L
	10 Liters	6001	2 x 10 L
	115 Liters	6003	Each
	210 Liters	6055	Each
Dri-Contrad® Powdered Autowasher Detergent	2 kG	2002	6 x 2 kG
	5 kG	2005	4 x 5 kG
	12.5 kG	2007	Each
	50 kG	2004	Each
Contrad® 100 Clean-In-Place Detergent	3.8 Liters	1504	4 x 3.8 L
	20 Liters	1503	Each
	115 Liters	1506	Each
	210 Liters	1505	Each



890 County Line Road
Bryn Mawr, PA 19010
(800) 332-6647
www.deconlabs.com