

Contrex™



Contrex™ is a comprehensive line of complimentary cleaning agents designed to meet any and all cleaning applications for labware, scientific apparatus and process equipment. (See product guide on reverse side).

Liquid

Powder

Concentrated

Alkaline

Low-foaming

Acidic

Cation-Free

Phosphate-Free

Enzymatic

Effervescent



460 Glennie Circle
King of Prussia, PA 19406
(800) 332-6647
www.deconlabs.com



Literature Numbers: BN0125056
BN0301051
DEC002

Contrex™ AP

Alkaline Powdered Detergent

Contrex AP is a concentrated powdered labware detergent formulated for use in manual, soak or ultrasonic cleaning of labware, scientific apparatus and process equipment. Contrex AP is completely soluble, free-rinsing and safe to dispose straight to drain. The user friendly bottle allows for easy dispensing and eliminates caking. Also available in a 1/2 ounce single serve which is pre-measured to make a 2 quart solution.

Contrex™ AL

Alkaline Liquid Detergent

Contrex AL is a concentrated alkaline liquid detergent formulated for use in manual, soak or ultrasonic cleaning of labware, scientific apparatus and process equipment. Contrex AL is free-rinsing, phosphate-free and disposable straight to drain.

Contrex™ AW

Powdered Autowasher Detergent

Contrex AW is a concentrated low-foaming powdered cleanser formulated for use in all automatic labware washing machines. Contrex AW is completely soluble, free-rinsing and biodegradable.

Contrex™ CA

Acidic Liquid Detergent

Contrex CA is a concentrated acidic liquid detergent made for use in descaling, mineral removal, oxide removal or in applications where an acidic detergent is indicated. Formulated with organic acid, Contrex CA can be used as a manual, soak or ultrasonic cleaner. It is phosphate-free and disposable straight to drain.

Contrex™ CF

Cation / Phosphate-Free Liquid Detergent

Contrex CF is a concentrated liquid detergent formulated for use in cleaning of equipment and electronics that are sensitive to metal cations and phosphates. Contrex CF is free-rinsing, phosphate-free and disposable straight to drain. Can be used for soak, scrub, sonication and auto washers.

Contrex™ EZ

Powdered Enzymatic Detergent

Contrex EZ is a powdered detergent formulated for use as an enzymatic detergent for use in manual or soak cleaning of labware and equipment contaminated with proteins and fats. Contrex EZ prevents corrosion, is free-rinsing and disposable straight to drain.

Contrex™ TL

Effervescent Detergent Tablets

Contrex TL Tablets are formulated for use in pipette washers, for parts cleaning or with any application assisted by effervescent action.

Contrex™ Product Guide

Item	Size	Decon Cat #	Fisher Cat #	Case Qty
Contrex™ AP <i>Alkaline Powdered Detergent</i>	50 x 1/2 oz. packets	5201	04-358-16	12pk/cs
	4 lbs.	5204	04-358-17	9
	25 lbs.	5225	04-358-18	—
	50 lbs.	5250	04-358-19	—
	100 lbs.	5200	04-358-20	—
Contrex™ AL <i>Alkaline Liquid Detergent</i>	1 qt.	5332	04-358-11	12
	1 gal.	5301	04-358-10	4
Contrex™ AW <i>Powdered Autowasher Detergent</i>	4 lbs.	5604	04-358-12	9
	25 lbs.	5625	04-358-13	—
	50 lbs.	5650	04-358-14	—
	100 lbs.	5600	04-358-15	—
Contrex™ CA <i>Acidic Liquid Detergent</i>	1 gal.	5701	04-358-21	4
	55 gal.	5755	04-358-22	—
Contrex™ CF <i>Cation / Phosphate-Free Liquid Detergent</i>	1 gal.	5801	04-358-26	4
Contrex™ EZ <i>Powdered Enzymatic Detergent</i>	4 lbs.	5404	04-358-23	9
	25 lbs.	5425	04-358-24	—
	50 lbs.	5450	04-358-25	—
Contrex™ TL <i>Effervescent Detergent Tablets</i>	100 tablets	5900	04-358-27	6