



# 200 Proof Ethanol

## 100 % Ethanol

### Technical Data Sheet

#### - BRIEF DESCRIPTION ó

Decon's 200 proof ethanol (100% ethanol) is a pure ethanol product made to USP specifications and used for research, healthcare, industrial, or production applications as a sanitizer, solvent, or reagent.

Please note that pure ethanol is regulated by the US ATF/TTB and depending on the taxable status of your institution, pure ethanol products are taxed according to the government schedule of taxes detailed in US CFR Part 27 (211.1119).

To purchase you must have an Industrial Alcohol User Permit from the AFT/TTB, specifically marked for the withdraw and use of alcohol free of tax. Tax free pure ethanol is intended for industrial use only. To apply for this permit visit [www.ttb.gov](http://www.ttb.gov) click Laboratory, click Laboratory Forms- click Form 5150.22

#### - INGREDIENTS-

100 % Ethyl alcohol CAS Number 64-17-5

#### - STORAGE -

Do not store Ethanol at temperatures above 120 degrees F.

#### - PRECAUTIONARY STATEMENT -

**Flammable** ó Keep away from heat, sparks, and flame. Surfaces that are hot may ignite even liquid product in the absence of sparks or flame.

**Caution** – May be harmful or fatal if swallowed. May irritate body tissues. Use with adequate ventilation. Avoid breathing vapor. Do not get in eyes, on skin, on clothing. Wash thoroughly after handling.

#### - PACKAGING -

- 1 Single Pint / each
- 1 Pint or 1 Glass Pint / 24 per case
- 1 Gallon or 1 Glass Gallon / 4 per case
- 5 Gallons / each
- 55 Gallons / each

#### - CATALOG NUMBERS -

<u>Manufacturer</u>	<u>Item</u>	<u>Cat. No.</u>
Decon Labs, Inc.	24x1 Pint	2716
	24x1 Pint, Glass	2716G
	1 Glass Single Pint	2716GEA
	4x1 Gallon	2701
	4x1 Gallon, Glass	2701G
	5 Gallon Plastic	2705
	5 Gallon Hedwin Cubitainer	2705HC
	5 Gallon Metal Container	2805M
	5 Gallon Spigot Closure	2705SG
	55 Gallon Drum	2755
	55 Gallon Metal Drum	2755M

For Additional Technical Information call  
Decon Customer Service  
(800) 332-6647, or  
See our website at [www.deconlabs.com](http://www.deconlabs.com)

Decon Labs, Inc. Revised 01/16/08