



# AQUAPUR™ ST

## Technical Data Sheet

---

### - BRIEF DESCRIPTION -

*AquaPur™ST is Purified Water that has been sterilized and packaged in a double bag for use in the cleanroom areas within pharmaceutical, biotechnology, and medical device facilities or anywhere the use of a sterile, USP Purified water is required.*

*AquaPur™ST is prepared and assayed to specifications which meet USP Purified Water Reagent Grade Water. AquaPur™ST meets CAP and NCCLS Standards for Reagent Grade Type 1 water at time of manufacture. AquaPur™ST contains no preservatives. It is gamma-irradiated to a sterility assurance level (SAL) of  $10^{-6}$ . Each lot is tested for sterility via USP 14-day sterility test. A copy of lot specific documentation, which details irradiation, sterility, and water quality limits for USP Purified Water, is attached to every case of AquaPur™ST.*

### - PRECAUTIONARY STATEMENT -

Not to be used for Parenteral Administration. Not for use in sterile irrigation, injection, infusion, or as a potable product. AquaPur™ ST is intended for manufacturing and general laboratory use.

### - STORAGE -

Store in cool area. Sterility can not be assured if the plastic bags are not sealed. AquaPur™ ST

has an expiration date of two years and the date is listed on the bottle.

### - PRODUCT SPECIFICATIONS -

AquaPur™ Sterile, USP grade Purified Water, is filtered through a 0.2  $\mu$ m filter. AquaPur™ ST is tested to meet USP specifications for Purified Water.

USP Water Conductivity  $\leq 1.3$  umhos/cm  
USP <645>

USP Total Organic Carbon < 0.5 PPM  
USP <643>

### - PACKAGING -

1 Gallon (3.8 Liters) / 4 per case

### - CATALOG NUMBER -

<u>Manufactured for:</u>	<u>Item</u>	<u>Cat. No.</u>
Decon Labs, Inc.	1 gallon	8904

---

**For Additional Technical Information call  
Decon Customer Service**

**(800) 332-6647, or**

**See our website at [www.deconlabs.com](http://www.deconlabs.com)**

---

®Decon Labs, Inc. Revised March 14, 2008  
(for verbage only, Feb 2008 version is also valid)