

- BRIEF DESCRIPTION -

FisherBrand's Bluing tablets are to be used to color the water in a toilet, via a blue dye, in order to prevent the water from being used to dilute a urine sample for drug testing. Packaged as two tablets per each foil packet, bluing tablets are an efficient and inexpensive means of preventing critical drug sample tampering.

- USE PROTOCOL -

Place one or two tablets in the toilet bowl. The effervescent action will immediately start to color the water blue.

- STORAGE -

Store the foil packets in a cool, dry place.

- PRECAUTIONARY STATEMENT -

DO NOT INGEST. For adding dye to toilet water only.

- PACKAGING -

2 tablets per Foil Pack/50 Packs per box

- CATALOG NUMBER -

<u>Manufacturer</u>	<u>Cat. No.</u>
FoxWest Sales	1220

**For Additional Technical Information call
FoxWest Sales Customer Service
(800) 332-6647, or
See our website at www.deconlabs.com**

Revised 11/16/06.