



CiDehol[®] ST Wipes

Technical Data Sheet

- BRIEF DESCRIPTION -

CiDehol ST Wipes are sterile wipers that are nonwoven polyester/cellulose wipes pre-saturated with 70% isopropyl alcohol for use in cleanrooms and sterile areas. Each pack is a convenient, re-sealable dispensing pouch containing 30 durable wipes, 9 x 9 inches. Each pack is double-bagged and gamma-irradiated to a sterility assurance level (SAL) of 10⁻⁶. A Lot Specific Document is shipped with each case. CiDehol ST wipes have been tested to show the number of particulates (equal to/greater than 0.5 um) from each wipe is less than 3.0 x 10⁴ per cm². The amount of extractables from water is less than 0.1 % and from isopropyl alcohol is less than 0.05 %.

- DIRECTIONS FOR DISPENSING -

The C-folded wipers are easily drawn from the re-sealable package and limits use to one wiper at a time.

- USE PROTOCOL -

Contents are sterile unless package has been damaged or opened. These wipes are designed for use in cleanroom areas within pharmaceutical, biotechnology and medical device facilities.

- ACTIVE INGREDIENTS -

Isopropyl Alcohol, USP Grade (CAS #67-63-0).....70 %
Deionized Water (CAS#7732-18-5).....30 %

- STORAGE AND DISPOSAL -

Store in cool area. Do not reuse wipes. Dispose of used wipe and empty dispenser in trash.

- PRECAUTIONARY STATEMENT -

CAUTION: FLAMMABLE. Do not expose wipes to heat, open flame, sparks, or elevated temperatures.

- FIRST AID -

For all cases of exposure consult with a Physician.
Eye Contact: Immediately flush with water for at least 15 minutes.
Skin Contact: Immediately rinse with water.
Inhalation: Immediately remove person to fresh air.
Ingestion: Immediately provide water to dilute ingested material.

- PACKAGING -

30 wipes per pouch / 24 pouches per case

- CATALOG NUMBER -

<u>Manufactured for:</u>	<u>Item</u>	<u>Cat. No.</u>
Decon Labs, Inc.	30 wipes	8730

**For Additional Technical Information call
Decon Customer Service
(800) 332-6647, or
See our website at www.deconlabs.com**

®Decon Labs, Inc. Revised 05/2003