



DECON'S PURE ETHANOL 200 PROOF



Decon's 200 proof ethanol (100% ethanol) is a pure ethanol product made to USP specifications and used for research, healthcare, industrial, or production applications as a cleaner, solvent or reagent.

It is available in convenient packaging configurations from stock. Custom solutions or packaging configuration may be available.

Please note that pure ethanol is regulated by the US ATF/TTB and depending on the taxable status of your institution, pure ethanol products are taxed according to the government schedule of taxes detailed in 27 CFR Chapter 1 section 27.41.

To purchase tax-free - you must have an Industrial Alcohol User Permit from the AFT/TTB, specifically marked for "the withdraw and use of alcohol free of tax" Tax free pure ethanol is intended for industrial use only. To apply for this permit visit www.ttb.gov and search for "Form 5150.22" - which is the application for the Industrial Alcohol User Permit.

STORAGE: Do not store Ethanol at temperatures above 120° F.

INGREDIENTS
100% Ethyl Alcohol (CAS #64-17-5)

Catalog Numbers / Packaging

Item	CASE QTY	Decon # Tax-Free	Decon # Tax-Paid
200 Proof, Pint (plastic)	24	2716	2716TP
200 Proof, Pint (glass)	24	2716G	2716GTP
200 Proof, 1 Single Pint (plastic)	-	2716GEA	2716GEATP
200 Proof, Gallon (plastic)	4	2701	2701TP
200 Proof, Gallon (glass)	4	2701G	2701GTP
200 Proof, 5 Gallon (plastic w/ spout)	-	2705	2705TP
200 Proof, 5 Gallon (Hedwin Cubitainer w/ spigot)	-	2705HC	2705HCTP
200 Proof, 5 Gallon (metal w/ spout)	-	2805M	2805MTP
200 Proof, 5 Gallon (w/ a spigot closure*)	-	2705SG	2705SGTP
200 Proof, 55 Gallon (plastic drum)	-	2755	2755TP
200 Proof, 55 Gallon (metal drum)	-	2755M	2755MTP

*Spigot available separately

If you possess an Industrial Alcohol Users Permit order using the Decon Tax-Free part number. All others would order using the Tax-Paid part number.

Rev. 08/2024