



INOCULATING NEEDLES



LeStab®, Decon's inoculating needles, are made of high quality, nichrome wire which is designed to cool rapidly after flaming. Most of our needles are permanently attached to lightweight, insulated handles which are perfectly balanced for easy manipulation. The nichrome wire cools rapidly after flaming.

Decon also offers 3 inch inoculating needles. These needles fit in Decon's LeLoop Brass holders. The Brass Holders come in 6 inch and 8 inch size.

USE PROTOCOL: LeStab®, Decon's inoculating needles with attached handles can be used and flamed many times before degrading the nichrome wire. Once degradation occurs, dispose of the needle with the attached handle and use a new one.

When the 3 inch inoculating needle needs to be replaced, remove the needle from the brass handle. Dispose of the used needle and place a fresh needle on the Brass holder.

Catalog Numbers / Packaging

Item	QTY	Decon #
MicroNeedle (11" w/ silver handle)	5/pack	MP186-5
Angled Needle (8" w/ blue handle)	1	MP4010-1G
MicroNeedle (8" w/ white handle)	1	MP4000-1G
Microneedles (3" - no handle)	12/pack	MP3325-12G
Brass Loop Holder 6"	10/pack	MP6-10
Brass Loop Holder 8"	10/pack	MP8-10

Rev. 12/2022