



# 200 Proof Molecular Grade Ethanol 100 % Ethanol

## Technical Data Sheet

---

### - BRIEF DESCRIPTION -

Decon's 200 proof ethanol (100% ethanol) is a Molecular Biology Grade, pure ethanol product made to USP specifications and is tested for DNAase / RNAase activity and its ability to precipitate DNA from solution. Molecular Biology Grade Ethanol can be used in absolute form or diluted for use in most applications in Molecular Biology Laboratories.

Pure Ethanol products are taxed according to the government schedule of excise taxes detailed in 27 CFR Chapter 1 section 27.41. This excise tax is included in the cost of this product. There is no need for Federal Permits.

### - INGREDIENTS-

100 % Ethyl alcohol    CAS Number 64-17-5

### - STORAGE -

Do not store Ethanol at temperatures above 120 degrees F.

---

### - PRECAUTIONARY STATEMENT -

**Flammable** – Keep away from heat, sparks, and flame. Surfaces that are hot may ignite even liquid product in the absence of sparks or flame.

**Caution** – May be harmful or fatal if swallowed. May irritate body tissues. Use with adequate ventilation. Avoid breathing vapor. Do not get in eyes, on skin, on clothing. Wash thoroughly after handling.

### - PACKAGING -

1 Pint Glass

### - CATALOG NUMBERS -

<u>Size</u>	<u>Cat. No.</u>
1 Pint	3916EA
1 Pint x 24 case	3916

---

**For Additional Technical Information call  
Decon Customer Service  
(800) 332-6647, or  
See our website at [www.deconlabs.com](http://www.deconlabs.com)**

---

Decon Labs, Inc.

Revised 5/18/2012