



# Platinum Loops Inoculating Loops Technical Data Sheet

## - BRIEF DESCRIPTION -

Decon's Platinum Loops are made of high quality, platinum wire which make these sturdy loops ideal for multiple re-use. The loops slide smoothly across the agar surface without tearing or gorging. Decon's Platinum loops are made to be used with Decon's brass holders.

## - USE PROTOCOL -

Decon's Platinum Loops can be used and flamed many times before degrading the platinum wire. Once degradation occurs, dispose of loop and attached a new loop to the handle.

## - APPLICATIONS -

Manufactured to critical tolerances, Decon's Platinum Loops can be used in quantitative procedures such as urine counts, etc.

## - PACKAGING -

Platinum Loop (Pt/15%Ir)2mm/.41/45	1 each	10/pk
Platinum Loop (Pt/15%Ir)3mm/.51/45	1 each	10/pk
Platinum Loop (Pt/3.5%RhPt/3.5%Rh)4mm/.91/25	1 each	10/pk
Platinum Loop (Pt/15%Ir)5mm/.51/45	1 each	10/pk

## - CATALOG NUMBERS -

<u>Manufacturer</u>	<u>Item</u>	<u>Cat. No.</u>
<b>Decon Labs</b>	Platinum Loop (Pt/15%Ir)2mm/.41/45 Fisher Cat # 13-102 A	MP202-1
	Platinum Loop (Pt/15%Ir)3mm/.51/45 Fisher Cat # 13-102 B	MP203-1
	Platinum Loop (Pt/3.5%RhPt/3.5%Rh)4mm/.91/25 Fisher Cat # 13-103	MP304-1
	Platinum Loop (Pt/15%Ir)5mm/.51/45 Fisher Cat # 13-102C	MP205-1

**For Additional Technical Information call  
Decon Customer Service Hotline  
(800) 332-6647, or  
See our website at [www.deconlabs.com](http://www.deconlabs.com)**

Decon Labs, Inc. *Revised 06/2006*